


* **BICO** *

- Same rugged performance as the UD, but belt drive.

- Motor equipped models with adjustable mounting blocks.

- Efficiency, simplicity of operation and maintenance.



- Samples ¼" or smaller to 150-200 mesh, depending on material.

- Adjusting screw with lock lever to hold permanent setting.

- Many models to choose from and various types of plates.

UA PULVERIZER V-BELT DRIVE

★ **242-53 UA V-BELT PULVERIZER.** 2 HP, 220 or 440 volt, 3 phase, 50 or 60 cycle motor.

★ **242-53X2S UA V-BELT PULVERIZER.**
2 HP, 110 or 220 volt, 1 phase, 50 or 60 cycle motor. MUST SPECIFY CYCLE.

★ **242-53X3 UA V-BELT PULVERIZER.** 3 HP, 220 or 440 volt, 3 phase, 50 or 60 cycle motor.

★ **242-54 UA V-BELT PULVERIZER WITH V-BELTS, MOTOR PULLEY, W/O MOTOR.**
Specify motor speed, shaft diameter, key-way size.

★ All UA Pulverizers supplied with 1 set iron alloy pulverizing plates ★

★ Safety guards available at extra charge ★

★ Special voltages and HP are available — Quoted on request ★

★ NOTE: thermal relays supplied for 220 volt operation unless otherwise specified ★



PO. BOX 6339
BURBANK, CA 91510

SINCE 1888

PHONE: 818-842-7179
FAX NO: 818-842-7976
TELEX: 194561 LSA-BICO